

USB Foot Switch

for Inspectis software, one pedal



Quality of Inspection

USB foot switch for hands free image capture in Inspectis software.



Foot switch is plugged into the USB port of computer and automatically recognized by the Inspectis software. Provides the user the freedom of using both hands when carrying out inspection tasks.

Specification:

Interface	USB 2.0
Cable length	1 800 mm
Colour	Black
Size	115x70x37 mm
Weight	210 g.

Ordering Number/Name

HD-051 USB Foot switch for Inspectis software, one pedal

Reseller:



info@microscopeworld.com | +1-760-438-0528

V1, 2018

www.inspect-is.com