

Aperture Window Plates

for resin mounted prepared samples



Quality of Inspection

Rigid metallic plates with stainless steel finish and 3 different aperture sizes for optical inspection of prepared samples with Inspectis Digital Inverted Microscopes. Compatible with all type of samples moulded in thermoplastic, epoxy, acrylic, conductive or other type of resins.



Glass window of the Inverted Microscope is simply replaced by any of the aperture window plates depending to size of the sample.

Specification:

Aperture sizes	Φ24mm, Φ38mm, Φ48mm
Material	Stainless steel plates. Bottom coated with black aluminium
Weight	100 g. max

Order Numbers:

HD-242 Inspectis Aperture Window Plates for resin mounted prepared samples (kit incl. 3 apertures)



MicroscopeWorld