



Cargille Laboratories
55 Commerce Rd - Cedar Grove, NJ 07009-1289
Phone : 973-239-6633 - Fax : 973-239-6096 - Web : www.cargille.com

Catalog # 1809

Typical Characteristics

Refractive Index Liquid **Series A**
Refractive Index **1.51000** at 589.3 nm and 25 ° C

Composition Aliphatic / Alicyclic Hydrocarbons and Hydrogenated Terphenyl
Appearance Colorless to light yellow liquid
Odor Slight
Color Stability In Direct Sun: may slightly darken after 1 day, more after 4 months, dark after 6 years
Index Change Rate Low : -0.00002 to +0.00005
by Evaporation expected after 32 days with exposed surface area to volume ratio of 0.2 cm²/cc @ 25 °C

Index at 20°C 1.51202
Pour Point < -7 °C
Boiling Point > 262 °C (760 mm Hg)
Flash Point > 138 °C (COC)
Brix Value (Per ICUMSA) 88.0 at 20°C
Density 0.910 g / cc at 25 °C
Density Temp Coef -0.0007 g / cc / °C
Coef of Thermal Expansion 0.0008 cc / cc / °C
Thermal Conductivity n/a cal / sec / cm² / °C - 1 cm thickness

Viscosity 35 cSt at 25 °C 32 cP at 25 °C
Surface Tension 31 dynes/cm at 25 °C

Soluble Carbon Tetrachloride, Ethyl Ether, Heptane, Methylene Chloride, Naphtha, Toluene, Turpentine, Xylene

Partly Soluble Acetone, Ethanol

Insoluble Water

Compatible 10 Month Immersion at 25 °C : Acrylic, Cellulose Acetate, Epoxy, Mylar, Nylon, Polycarbonate, Polyester, Polyethylene, Polypropylene, Polyurethane, Polyvinyl Chloride, Phenolic, Teflon, Silicone (Sylgard 184) and Fluorosilicone (730 RTV) Rubbers, Tygon F-4040-A, Tygothane; Aluminum, Brass, Copper, Steel

Incompatible Polystyrene, Latex, Neoprene, Silicone (3140 RTV) Rubber, Tygon S-50-HL, R-3603, B-44-3



Refractive Index Liquid **Series A**
 Refractive Index **1.51000** at 589.3 nm and 25 °C

Cauchy Coefficients
 A 1.492203E+00 B 5.7180E+03 C 1.6058892E+08

Cauchy Equation at 25°C $A + B / \lambda^2 + C / \lambda^4$ (λ = Wavelength in nm)

| Wavelength (nm) | Refractive Index | Transmittance | | | |
|-----------------|------------------|---------------|------|------|-------|
| | | 0.1 mm | 1 mm | 1 cm | 10 cm |
| 225.0 | - | - | - | - | - |
| 240.0 | - | - | - | - | - |
| 250.0 | - | - | - | - | - |
| 270.0 | - | - | - | - | - |
| 290.0 | - | - | - | - | - |
| 308.0 | - | - | - | - | - |
| 355.0 | 1.548 | 86 | 23 | 0 | 0 |
| 365.0 | 1.544 | 98 | 82 | 14 | 0 |
| 406.0 | 1.5328 | 100 | 95 | 63 | 1 |
| 473.0 | 1.5210 | 100 | 99 | 94 | 56 |
| 532.0 | 1.5144 | 100 | 100 | 98 | 81 |
| 589.3 | 1.5100 | 100 | 100 | 99 | 89 |
| 632.8 | 1.5075 | 100 | 100 | 99 | 93 |
| 656.3 | 1.5063 | 100 | 100 | 100 | 95 |
| 790.0 | 1.5018 | 100 | 100 | 100 | 99 |
| 828.0 | 1.5009 | 100 | 100 | 100 | 97 |
| 981.0 | 1.498 | 100 | 100 | 96 | 66 |
| 1310.0 | 1.496 | 100 | 99 | 89 | 32 |
| 1550.0 | 1.495 | 100 | 98 | 83 | 15 |
| 2500.0 | - | - | - | - | - |
| 3700.0 | - | - | - | - | - |

($n_F - n_C$) 0.0129

Abbe v_D 39.4

Temp. Coefficient -0.000404 dn_D/dt (15 - 35 °C)

Shelf Life: 5 Years from Date of Manufacture for Unopened Bottles, Half the Remaining Time after Opening