



**Cargille Laboratories**  
**55 Commerce Rd - Cedar Grove, NJ 07009-1289**  
**Phone : 973-239-6633 - Fax : 973-239-6096 - Web : www.cargille.com**

Catalog # 1803

Typical Characteristics

Refractive Index Liquid      **Series AAA**  
Refractive Index              **1.30000**      at    589.3 nm    and      25      ° C

---

Composition                      Perfluorocarbon and Chlorofluorocarbon ( not the types thought to affect the ozone )  
Appearance                      Colorless Liquid  
Odor                                      None  
Color Stability                      In Direct Sun: no visible change after 10 years  
Index Change Rate                Moderate : -0.00030 to +0.00009  
by Evaporation                    expected after 32 days with exposed surface area to volume ratio of 0.2 cm<sup>2</sup>/cc @ 25 °C

---

Index at 20°C                      1.30167  
Pour Point                              < -20 °C  
Boiling Point                        > 215 °C ( 760 mm Hg )  
Flash Point                            None °C ( COC )  
Brix Value (Per ICUMSA)            NA  
Density                                1.939 g / cc at 25 °C  
Density Temp Coef                  -0.0022 g / cc / °C  
Coef of Thermal Expansion        0.0011 cc / cc / °C  
Thermal Conductivity                0.00017 cal / sec / cm<sup>2</sup> / °C - 1 cm thickness  
  
Viscosity                              14 cSt at 25 °C                              27 cP at 25 °C  
Surface Tension                      18 dynes/cm at 25 °C

---

Soluble                                Other Chlorofluorocarbons, Galden PFS2, Fluoroclean HE

Partly Soluble                      Most organic solvents

Insoluble                              Water

---

Compatible                        10 Month Immersion at 25 °C: Acrylic, Cellulose Acetate, Epoxy, Mylar, Nylon, Polycarbonate, Polyester, Polyethylene, Polypropylene, Polystyrene, Polyurethane, Polyvinyl Chloride, Phenolic, Teflon, Latex, Neoprene, Silicone ( Sylgard 184, 3140 RTV ), and Fluorosilicone ( Silastic 730 RTV ) Rubbers ; Tygon F-4040-A, Tygothane, Brass, Copper, Steel

Incompatible                      Burna-S, Natural, and some Silicone Rubbers; Tygon S-50-HL, R-3603, B-44-3; Chlorotrifluoro Ethylene Polymers, Aluminum

---



Refractive Index Liquid **Series AAA**  
 Refractive Index **1.30000** at 589.3 nm and 25 °C

Cauchy Coefficients  
 A 1.293552E+00 B 2.6171E+03 C -1.3126839E+08

Cauchy Equation at 25°C  $A + B / \lambda^2 + C / \lambda^4$  ( $\lambda$  = Wavelength in nm)

Wavelength (nm)	Refractive Index	Transmittance			
		0.1 mm	1 mm	1 cm	10 cm
225.0	-	-	-	-	-
240.0	1.30	99	94	52	0
250.0	1.30	100	96	66	1
270.0	1.30	100	99	91	37
290.0	1.31	100	100	96	64
308.0	1.31	100	100	98	80
355.0	1.306	100	100	99	92
365.0	1.306	100	100	99	94
406.0	1.3046	100	100	100	96
473.0	1.3026	100	100	100	99
532.0	1.3012	100	100	100	99
589.3	1.3000	100	100	100	99
632.8	1.2993	100	100	100	99
656.3	1.2989	100	100	100	99
790.0	1.2974	100	100	100	99
828.0	1.2971	100	100	100	95
981.0	1.296	100	100	100	95
1310.0	1.295	100	100	100	95
1550.0	1.295	100	100	98	82
2500.0	1.29	100	99	87	25
3700.0	1.29	97	70	3	0

( $n_F - n_C$ ) 0.0034

Abbe  $v_D$  89.4

Temp. Coefficient -0.000333  $dn_D/dt$  (15 - 35 °C)

Shelf Life: 5 Years from Date of Manufacture for Unopened Bottles, Half the Remaining Time after Opening