

## REFRACTIVE INDEX LIQUID SERIES E

n( 5893 Å) 25°C =1.6150

## TYPICAL CHARACTERISTICS

COMPOSITION ..... Hydrogenated Terphenyl and  
1-Bromonaphthalene

APPEARANCE ..... Light yellow liquid

ODOR ..... Slight: unpleasant

COLOR STABILITY ..... In sun: may slightly darken after 1 day,  
becoming very dark after 4 months, and dark with precipitate after 6 years

INDEX CHANGE RATE BY EVAPORATION . Low: -0.00003 to +0.00003 expected:  
exposed surface area to volume ratio of 0.2 cm<sup>2</sup>/cc @ 25°C for 32 days

POUR POINT °C ..... <6

BOILING POINT °C @ 760mm Hg ..... >279

FLASH POINT °C CC ..... >113

DENSITY g/cc @ 25°C ..... 1.253

DENSITY TEMP. COEF. g/cc/°C ..... -0.0008

COEF. OF THERM. EXP. cc/cc/°C .... 0.0007

VISCOSITY centistokes @ 25°C ..... 18 (ca. 29 @ 15°C, 14 @ 35°C)

SURFACE TENSION dynes/cm @ 25°C .. 39

SOLUBLE: Acetone, Carbon Tetrachloride, Ethyl Ether, Freon TF, Heptane,  
Methylene Chloride, Naphtha, Toluene, Turpentine, Xylene

PARTLY SOLUBLE: Ethanol; INSOLUBLE: Water

COMPATIBLE 9 month immersion @ 25°C: Acrylic, Cellulose Acetate, Epoxy,  
Mylar, Nylon, Polyester, Polyethylene, Polypropylene, Polyurethane,  
Polyvinyl Chloride, Phenolic, Teflon; Silicone (Sylgard 184, 3140 RTV)  
and Fluorosilicone (Silastic 730 RTV) Rubbers; Tygothane; Aluminum, Steel  
(tests done on one example of each)

INCOMPATIBLE: Polycarbonate, Polystyrene, Latex, Neoprene, Tygon (all types  
except Tygothane), (Acrylic @ 55°C). May tarnish Copper and Brass

TOXICITY ..... Moderate in our experience (request MSDS)

CAUCHY EQUATION: refractive index as a function of wavelength at 25°C

W = wavelength in angstroms (Å)

$$n(W) = 1.582834 + ( 935849.7 )/W^2 + ( 6.29917E+12 )/W^4$$

SOURCE OR SPECTRAL LINE	WAVELENGTH (angstroms)	REFRACTIVE INDEX 25°C	% TRANSMITTANCE 25°C		
			1mm	1cm	10cm
near UV cut off	3500	1.701	40	0	0
i (Hg)	3650	1.689	76	6	0
h (Hg)	4047	1.6635	93	47	0
F' (Cd)	4800	1.6353	99	87	25
F (H)	4861	1.6337	99	88	28
e (Hg)	5461	1.6213	100	97	70
D (Na D1,D2 mean)	5893	1.6150	100	98	82
HeNe laser	6328	1.6101	100	98	86
C' (Cd)	6439	1.6091	100	98	86
C (H)	6563	1.6080	100	99	88
Ruby laser	6943	1.6050	100	99	89
GaAs laser	8400	1.5974	100	99	92
Nd:YAG laser	10648	1.592	100	98	84
Diode	13000	1.589	99	93	47
Diode	15500	1.587	99	89	31

$n_F - n_C$  = 0.0258

Abbe  $v_D$ :  $(n_D - 1)/(n_F - n_C)$  = 23.9

Temp. coef:  $dn_D/dt$  15-35°C = -0.000446

CARGILLE LABORATORIES INC.

55 COMMERCE RD. CEDAR GROVE N.J. 07009-1289 U.S.A.

phone: (973)239-6633

Fax: (973)239-6096